


TYPICAL DRAWINGS
VPL UNENCLOSED NO PIT
ADJ/90° EXIT PLATFORM



PRODUCT SPECIFICATIONS

RATED LOAD:	750 LBS (340 KG) MAXIMUM
INPUT POWER SOURCE:	
-AC POWERED UNIT:	110-120 VAC, 20 AMP, 60 HZ DEDICATED SERVICE
-OPT BATTERY BACKUP UNIT:	BASE 2X 12 V BATTERIES (NOT AVAILABLE FOR VPL800-1400) PLUS 4X 12 V BATTERIES
DRIVE:	
-AC POWERED UNIT:	2/3 HP, 26:1 GEARMOTOR W/BRAKE, 90 VDC, INTERMITTENT DUTY, CLASS H
DRIVE:	1" DIAMETER, 0.1" PITCH, 10 START ACME SCREW (1" LEAD)
LIFT SPEED:	7-10 FEET/MINUTE
SAFETY:	MOTOR FRICTION BRAKE AND ELECTRIC BRAKE, OVERSPEED GOVERNOR, BACKUP ACME NUT, SAFETY PAN, E-STOP, OVERLOAD SENSING, OVER TEMPERATURE SENSING, PERFORMANCE MONITORING, BATTERY MONITORING (ON BATTERY BACKUP UNITS)

SPECIFICATIONS SUBJECT TO CHANGE WITH OR WITHOUT NOTICE.

DESIGN AND REGULATORY REQUIREMENTS

USA FOOD & DRUG ADMINISTRATION: CLASS II, 510(K) EXEMPT FILE NO. 890.3930, PRODUCT CODE: PCE

- DESIGNED IN COMPLIANCE WITH THE FOLLOWING CODES:
- ASME A18.1 SAFETY STANDARD FOR PLATFORM LIFTS AND STAIRWAY CHAIRLIFTS (SECTIONS 2 AND 5)
 - CSA B44.1/ASME A17.5 ELEVATOR AND ESCALATOR ELECTRICAL EQUIPMENT
 - ICC ANSI A117.1 ACCESSIBLE AND USABLE BUILDINGS AND FACILITIES
 - NFPA 70 NATIONAL FIRE PROTECTION AGENCY - NATIONAL ELECTRICAL CODE
 - ETL CERTIFIED - CONTROL NUMBER 3148125


FDA INDICATION OF USE STATEMENT

THE HIGHLANDER II VERTICAL PLATFORM LIFT IS TO ASSIST IN THE TRANSFER OF PATIENTS OR A MOBILITY IMPAIRED PERSON, UP AND DOWN BETWEEN LEVELS OF A RESIDENTIAL OR COMMERCIAL FACILITY.



CONFIDENTIAL INFORMATION NOT TO BE COPIED OR USED WITHOUT WRITTEN CONSENT. INTERPRET DRAWING AND TOL PER ASME Y14.5 - 2009		FREEOM Lift Systems	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES ARE:		DWG NO.	REV
DECIMALS .X = ±0.1 .XX = ±0.02 .XXX = ±0.010		ENG-000842	A
HOLE DIA'S = ±.005 ANGLES = ±1°		DESCRIPTION TYPICAL DRAWINGS VPL UNENCLOSED NO PIT ADJ/90° EXIT PLATFORM	
DO NOT MEASURE FROM DRAWING		DRAWN BY C POLSTER	
		DATE 17FEB20	
		SCALE NTS	SHEET 1 OF 4

TYPICAL DRAWINGS
VPL UNENCLOSED NO PIT
ADJ/90° EXIT PLATFORM

REVISIONS				
REV.	DESCRIPTION	BY	CHK	DATE
-	See Sheet1	-	-	-



RIGHT HAND UNIT SHOWN



NOTES:

- 1 BOTH RIGHT AND LEFT HAND PLATFORMS AVAILABLE
- 2 PLATFORM GATE ADDS 1.50" TO DIMENSIONS E & H
- 3 ASME A18.1 REQUIRES COMMERCIAL LIFTS OVER 60" TRAVEL TO BE IN AN ENCLOSED HOISTWAY.

PLATFORM DIMENSIONS TABLE							
PLATFORM	D	E	F	G	H	J	K
36" X 54"	36.00	54.00	9.60	33.44	78.4	33.33	54.78
36" X 60"	36.00	62.00	13.60	33.44	84.4	33.33	54.78
42" X 54"	42.00	54.00	9.60	39.44	78.4	36.33	60.78
42" X 60"	42.00	62.00	13.60	39.44	84.4	36.33	60.78

TOWER DIMENSIONS TABLE			
LIFT MODEL	A	B	C
VPL400	77	53	95
VPL600	101	77	119

CONFIDENTIAL INFORMATION
NOT TO BE COPIED OR USED
WITHOUT WRITTEN CONSENT.
INTERPRET DRAWING AND TOL
PER ASME Y14.5 - 2009

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES.
TOLERANCES ARE:
DECIMALS .X = ±0.1
.XX = ±0.02
.XXX = ±0.010
HOLE DIA'S = ±.005
ANGLES = ±1°
DO NOT MEASURE FROM DRAWING

FREEOM Lift Systems

DWG NO. **ENG-000842** REV **A**


DESCRIPTION
**TYPICAL DRAWINGS VPL UNENCLOSED
NO PIT ADJ/90° EXIT PLATFORM**

SCALE NTS SIZE: B SHEET 2 OF 4


DRAWN BY **C POLSTER** DATE **17FEB20**

TYPICAL DRAWINGS
VPL UNENCLOSED NO PIT
ADJ/90° EXIT PLATFORM


REVISIONS				
REV.	DESCRIPTION	BY	CHK	DATE
-	See Sheet1	-	-	-




PAD DETAILS





LIFT WITH AUTO RAMP



LIFT WITH GATE AND AUTO RAMP



SECTION VIEW
OF CONCRETE ANCHORS
(NOT TO SCALE)



PAD DIMENSIONS MINIMUM					
PLATFORM	N	P	R	T	U
36" X 54"	8.44	56.18	78.4	79.9	93.82
36" X 60"	12.44	56.18	86.4	87.9	101.82
42" X 54"	8.44	62.18	78.4	79.9	99.82
42" X 60"	12.44	62.18	86.4	87.9	107.82

CONFIDENTIAL INFORMATION
NOT TO BE COPIED OR USED
WITHOUT WRITTEN CONSENT.
INTERPRET DRAWING AND TOL
PER ASME Y14.5 - 2009

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES.
TOLERANCES ARE:

DECIMALS .X = ±0.1
 .XX = ±0.02
 .XXX = ±0.010

HOLE DIA'S = ±.005

ANGLES = ±1°

DO NOT MEASURE FROM DRAWING

DWG NO.
ENG-000842

REV
A

DESCRIPTION
**TYPICAL DRAWINGS VPL UNENCLOSED
NO PIT ADJ/90° EXIT PLATFORM**

DRAWN BY
C POLSTER


DATE
17FEB20

SCALE NTS SIZE: B

SHEET 3 OF 4

TYPICAL DRAWINGS
VPL UNENCLOSED NO PIT
ADJ/90° EXIT PLATFORM

REVISIONS				
REV.	DESCRIPTION	BY	CHK	DATE
-	See Sheet1	-	-	-



Site Construction Details

Electrical Requirements:

Check NFPA 70 and all applicable codes for all electrical and wiring requirements.

Platform Pathway Requirements:

Ensure the platform pathway is sufficiently illuminated, clear of any electrical conduit and wireways and there is sufficient headroom clearance (minimum of 79"-2007mm) per ASME A18.1 throughout floor-to-floor travel.

Floor Recommendations:

A 4" (102mm) thick, 3500 PSI minimum compressive strength, reinforced concrete pad is recommended. Refer to sheet 3 for minimum pad dimensions.

Tower Floor Anchoring:

VPL must be fastened to concrete pad using a minimum of four (4) 3/8" (5/16" bolt) x minimum 2-1/2" long concrete anchors suitable for the environment. Refer to Pad Details on sheet 3 for mounting hole locations. Follow selected concrete anchor manufacturer's guidelines and all applicable codes.

Tower Support:

Tower support is not required on lifts with travel under 8 feet. Lifts with travel of 8 feet or more are required to be supported appropriately. See views at left or installation manual for more details.

Top Landing Gate Attachment:

Refer to landing gate detail pages.

Platform to Top of Landing Sill Clearance:

ASME A18.1code stipulates the platform-to-sill clearance at the upper landing shall not be less than 3/8" (9.5mm) nor exceed 3/4" (19.1mm). Follow all applicable codes.


Fascia Wall Requirements:

ASME A18.1code stipulates that fascia should be smooth and/or non-perforated that guards the full length and width of the platform. The fascia shall be securely fastened from the upper landing sill down to ground level. The fascia must be able to support a 125-pound side load over any 4-inch square area. Follow all applicable codes.

CONFIDENTIAL INFORMATION NOT TO BE COPIED OR USED WITHOUT WRITTEN CONSENT. INTERPRET DRAWING AND TOL PER ASME Y14.5 - 2009		FREEOM Lift Systems	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES ARE: DECIMALS .X = ±0.1 .XX = ±0.02 .XXX = ±0.010 HOLE DIA'S = ±.005 ANGLES = ±1° DO NOT MEASURE FROM DRAWING		DWG NO. ENG-000842	REV A
DESCRIPTION TYPICAL DRAWINGS VPL UNENCLOSED NO PIT ADJ/90° EXIT PLATFORM		DRAWN BY C POLSTER	DATE 17FEB20
SCALE NTS		SIZE: B	SHEET 4 OF 4

TYPICAL DRAWING
LANDING GATES

REVISIONS				
REV.	DESCRIPTION	BY	CHK	DATE
A	INITIAL RELEASE, ECO-005013	CJP	DN	16APR20



- NOTES:
- 1 LEFT HAND LANDING GATE SHOWN, RIGHT HAND LANDING GATE MIRRORED. LANDING GATE IS FIELD CONVERTIBLE, IF INTERLOCK IS REPLACED.
 - 2 APPROPRIATE LANDING GATE SUPPORT FRAMING REQUIRED

LANDING GATE DIMENSIONS TABLE						
GATE WIDTH	A	B	C	D	E	F
36" WIDE	36.13	42.99	36.12	41.62	19.81	18.00
42" WIDE	42.13	48.99	42.12	47.62	22.81	21.00

CONFIDENTIAL INFORMATION
NOT TO BE COPIED OR USED
WITHOUT WRITTEN CONSENT.
INTERPRET DRAWING AND TOL
PER ASME Y14.5 - 2009

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES.
TOLERANCES ARE:

DECIMALS .X = \pm .1
 .XX = \pm 0.02
 .XXX = \pm 0.010

ANGLES = \pm 1°
HOLES = \pm .005

DO NOT MEASURE FROM DRAWING

FREEOM Lift Systems

DWG NO. **ENG-000847** REV **A**

DESCRIPTION

**TYPICAL DRAWINGS
LANDING GATES**

DRAWN BY: **C POLSTER** DATE: **17FEB20**

SCALE: **NTS** SIZE: **B** SHEET **1 OF 1**



Freedom Commercial Wheelchair Platform Lift Purchase Terms

To our knowledge the lift as outlined on quote # _____, drawing # _____ will meet ASME A18.1 2014 code when properly installed as per the code.

Purchaser's Responsibilities:

- 1) Freedom Commercial Wheelchair Platform Lifts are designed and configured to be ASME A18.1 compliant however the purchaser is responsible for confirming compliance with all other federal, state and local jurisdictions codes and/or regulations.
- 2) The lift must be installed by a qualified installer. Several states require a licensed elevator mechanic or technician to install the lift. The purchaser must confirm installation requirements for their state and/or local jurisdiction. If a licensed elevator mechanic or technician is not required, a licensed contractor is required to install the lift and a licensed electrician is required to hardwire the electrical.
- 3) Proof of liability insurance is necessary to purchase a commercial lift. At the time of purchase the purchaser must provide the installer's current certificate of insurance, evidencing a minimum of \$1,000,000 Commercial General Liability.

Accessibility Professionals Inc. will not be responsible for any costs associated with modifying the Commercial Wheelchair Platform Lift if deemed not complaint upon inspection by your local inspector and/or permitting office. In addition, please be advised that we are unable to accept returns based on orders placed incorrectly. For complete list of our terms and conditions please visit our web site.

<https://freedomshowers.com/terms-conditions>

In order to proceed with the Commercial Wheelchair Platform Lift order, we require a witnessed signature below to confirm your understanding of the requirements and responsibilities as the purchaser of the lift.

I _____ agree that I have read and understood the requirements and responsibilities as the purchaser of the Commercial Wheelchair Platform Lift.

Print Name Signature

This _____ day of _____, Year _____

I _____ WITNESS that the above purchaser has read and signed this form.

Print Name Signature

This _____ day of _____, Year _____